
SANOSIL
DISINFECTANTS FOR LIFE 

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Bremen ,den 28.01.2014

Summary: virus-inactivating properties of Sanosil S010 of Sanosil AG according to EN 14476:2013

This summary is based on the following test reports of MikroLab GmbH for the surface disinfectant Sanosil S010 produced by Sanosil AG for surface disinfection:

poliovirus test report 27.01.2014
adenovirus test report 27.01.2014
MNV test report 28.01.2014

The following concentration and exposure time are necessary for the inactivation of the 3 test viruses:

undiluted 30 minutes

in order to achieve a four log₁₀ reduction (inactivation ≥ 99.99 %) under clean conditions in a quantitative suspension test according to EN 14476:2013.

After evaluation with poliovirus type 1, adenovirus type 5 and murine norovirus (MNV) the surface disinfectant Sanosil S010 can be declared as having “virucidal activity” properties according to EN 14476:2013.

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of the product:

Sanosil S010

to be intended for: **Disinfection of surfaces used in the medical area**

The testing of the product was carried out according to

DIN EN 13727


Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of bactericidal activity in the medical area - Test method and requirements (phase 2, step 1); German version EN 13727:2012+A2:2015

This expert opinion is based on the test report "Prüfbericht PL 17-72 EN 13727 170807", dating 07 August 2017.

Results

The product meets the standards given by the DIN EN 13727 for the bactericidal efficacy under high organic burden and the following relations of concentration and time of action:

15 min time of action and 80 % concentration.

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of the product:

Sanosil S010

to be intended for: **Disinfection of surfaces used in the medical area**

The testing of the product was carried out according to

DIN EN 13624


Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of fungicidal or yeasticidal activity in the medical area - Test method and requirements (phase 2, step 1); German version EN 13624:2013

This expert opinion is based on the test report "Prüfbericht PL 17-72 EN 13624 170807", dating 07 August 2017.

Results

The product meets the standards given by the DIN EN 13624 for the yeasticidal efficacy under dirty conditions using the following relations of concentration and time of action:

**15 min time of action and 97 % concentration as well as
30 min time of action and 80 % concentration.**

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*EPA registered claims summary at time of printing. Subject to change.

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Comparison of disinfectants for surfaces in healthcare

Rating	★★★★★ Sodium hypochlorite	★★★★★ Alcohols	★★★★★ QACs	★★★★★ Aldehydes	★★★★★ Peracetic acid	★★★★★ Glutaraldehyde	★★★★★ Guanidines	★★★★★ Sanosil 5003	★★★★★ Sanosil 5010
Degradability/environmental compatibility	MODERATE (AOX)	GOOD	POOR	POOR	GOOD	POOR	POOR	GOOD	GOOD
Depositing effect	NO	NO	NO	NO	NO	NO	NO	YES	YES
Enhances antibiotic resistance	NO	NO	YES	NO	NO	NO	NO	NO	NO
Corrosion problems	YES	NO	NO	YES	YES	NO	NO	RARELY	RARELY
Resistance development	NO	NO	YES	YES	NO	YES	YES	NO	NO
Contaminant capacity	MODERATE	GOOD	MODERATE	MODERATE	MODERATE	MODERATE	MODERATE	MODERATE	MODERATE
Dermatological compatibility	POOR	GOOD	GOOD	POOR	POOR	POOR	POOR	GOOD	MEDIUM
Application on large surfaces	YES	NO	YES	YES	YES	YES	YES	YES	YES
Bactericidal + yeasticidal	YES	YES	YES	YES	YES	YES	YES	YES	YES
Fungicidal	MODERATE	MODERATE	GOOD	MODERATE	GOOD	MODERATE	MODERATE	MODERATE	GOOD
Sporicidal	YES	NO	YES	NO	YES	YES	NO	NO	YES
Virucidal	YES	PARTIALLY	YES	YES	YES	NO	NO	PARTIALLY	YES
Score comparison	15	16	15	10	17	12	10	18	22

Score: +2 +1 0


Comparison of disinfectants for surfaces in agriculture

Rating	★★★★★ Sodium hypochlorite	★★★★★ Phenols	★★★★★ QACs	★★★★★ Aldehydes	★★★★★ Peracetic acid	★★★★★ Glutaraldehyde	★★★★★ Guanidines	★★★★★ Sanosil 5010
Degradability/environmental compatibility	MODERATE (AOX)	POOR	MODERATE	POOR	GOOD	POOR	POOR	GOOD
Depositing effect	NO	NO	NO	NO	NO	NO	NO	YES
Enhances antibiotic resistance	NO	NO	YES	NO	NO	NO	NO	NO
Corrosion problems	YES	YES	NO	YES	YES	NO	NO	RARELY
Resistance development	NO	YES	YES	YES	NO	YES	YES	NO
Contaminant capacity	MODERATE	GOOD	MODERATE	MODERATE	MODERATE	MODERATE	MODERATE	MODERATE
Dermatological compatibility	POOR	POOR	MODERATE	POOR	POOR	POOR	POOR	MODERATE
Bactericidal + yeasticidal	YES	YES	YES	YES	YES	YES	YES	YES
Fungicidal	MODERATE	MODERATE	GOOD	MODERATE	GOOD	MODERATE	MODERATE	GOOD
Sporicidal	YES	PARTIALLY	YES	NO	YES	YES	NO	YES
Virucidal	YES	NO	YES	YES	YES	NO	NO	YES
Against oocysts and worm eggs	NO	YES	NO	NO	NO	NO	NO	NO
Score comparison	13	10	13	8	15	10	8	20

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CH-3003 Bern, FOPH **A-Priority**

Sanosil AG
Eichtalstrasse
CH- 8634 Hombrechtikon

Reference:
Your Ref.:
Our Ref.: BLO / sik
Liebefeld, 28th March 2018

Federal Department of Home Affairs FDHA
Federal Office of Public Health FOPH
Notification authority for chemicals

Certificate of Registration in Switzerland

We confirm that the following biocidal product is authorised in Switzerland

Name of product	SZID Number	Authorisation Number
Sanosil Super 25	150868	CHZB2260

This biocidal product is allowed to be placed on the Swiss market by Sanosil AG, pursuant to article 7 paragraph 1 letter d and article 13 of the Swiss Ordinance on Biocidal Products (OBP; SR 813.12).

Kind regards


Notification authority for chemicals
Scientific assistant

Dr. Olivier Blaser

Federal Office of Public Health
Postal address: CH-3003 Bern
www.bag.admin.ch, www.cheminfo.ch

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CH-3003 Bern, FOPH **A-Priority**

Sanosil AG
Eichtalstrasse
CH- 8634 Hombrechtikon

Reference:
Your Ref.:
Our Ref.: BLO / sik
Liebefeld, 28th March 2018

Federal Department of Home Affairs FDHA
Federal Office of Public Health FOPH
Notification authority for chemicals

Intended Usage

Product types

- 02-09 Surface disinfection general
- 02-12 Surface disinfection: Air conditioning, refrigeration, ventilation
- 02-26 Treatment of waste water, waste, chemical toilet
- 03-01 Animal husbandry: Preventive disinfection (define the species in tab "Comments")
- 02-11 Industrial surface disinfection
- 02-18 Disinfectants: Anti-mould
- 04-01 Disinfectants for food industry (incl. kitchens, restaurants, canteens): Surfaces
- 04-04 Disinfectants for beverage industry and brewery
- 05-01 Drinking water disinfection
- 02-10 Disinfectants: Surfaces for medical use, schools, pharmaceuticals, etc.

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Outbreak of a new Coronavirus in China

In the central Chinese city of Wuhan, there has been an extraordinary accumulation of pneumonia. These are caused by a previously unknown coronavirus (2019-nCoV).

Coronaviridae belong to the enveloped viruses and are therefore not particularly resistant to disinfectants.

Sanosil Disinfectants S003, SanoClean AR, S006, S010 and S015 are effective against Coronavirus 2019-nCoV.

The Coronaviridae are a virus family within the order Nidovirales. Their representatives cause very different diseases in different vertebrates such as mammals, birds and fish. Coronaviruses are genetically highly variable, and individual virus species can infect several host species by crossing the species barrier.

Such species transitions have led to infections in humans with, among other things, the SARS-associated coronavirus (SARS-CoV) – the pathogen that caused the 2002/2003 SARS pandemic – and the Middle East respiratory syndrome coronavirus (MERS-CoV), which appeared in 2012.

A cluster of pneumonia observed in the Chinese city of Wuhan around the turn of the year 2019/2020 is also attributed to a previously unknown coronavirus, which was given the provisional name 2019-nCoV. In humans, the human coronavirus species is of particular importance as a pathogen causing mild respiratory infections (colds) and severe acute respiratory syndrome.

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Persistence of coronaviruses on inanimate surfaces and their inactivation with biocidal agents

G. Kampf^{a,*}, D. Todt^b, S. Pfaender^b, E. Steinmann^b

^a University Medicine Greifswald, Institute for Hygiene and Environmental Medicine, Ferdinand-Sauerbruch-Straße, 17475 Greifswald, Germany

^b Department of Molecular and Medical Virology, Ruhr University Bochum, Universitätsstrasse 50, 44801 Bochum, Germany

ARTICLE INFO

Article history:
Received 31 January 2020
Accepted 31 January 2020
Available online xxx

Keywords:
Coronavirus
Persistence
Inanimate surfaces
Chemical inactivation
Biocidal agents
Disinfection

SUMMARY

Currently, the emergence of a novel human coronavirus, SARS-CoV-2, has become a global health concern causing severe respiratory tract infections in humans. Human-to-human transmissions have been described with incubation times between 2–10 days, facilitating its spread via droplets, contaminated hands or surfaces. We therefore reviewed the literature on all available information about the persistence of human and veterinary coronaviruses on inanimate surfaces as well as inactivation strategies with biocidal agents used for chemical disinfection, e.g. in healthcare facilities. The analysis of 22 studies reveals that human coronaviruses such as Severe Acute Respiratory Syndrome (SARS) coronavirus, Middle East Respiratory Syndrome (MERS) coronavirus or endemic human coronaviruses (HCoV) can persist on inanimate surfaces like metal, glass or plastic for up to 9 days, but can be efficiently inactivated by surface disinfection procedures with 62–71% ethanol, 0.5% hydrogen peroxide or 0.1% sodium hypochlorite within 1 minute. Other biocidal agents such as 0.05–0.2% benzalkonium chloride or 0.02% chlorhexidine digluconate are less effective. As no specific therapies are available for SARS-CoV-2, early containment and prevention of further spread will be crucial to stop the ongoing outbreak and to control this novel infectious threat.

© 2020 The Healthcare Infection Society. Published by Elsevier Ltd. All rights reserved.

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Swiss Confederation

Federal Department of Home Affairs FDHA

Institute of Virology and Immunology IVI
Biosafety and Engineering

CH-3147 Mittelhäusern, IVI

Sanosil AG
Herr D. Rüegg
Eichtalstrasse 49
8634 Hombrechtikon
Switzerland

Your reference:
Reference no. /File no.:
Our reference: KUD
Contact person: Daniel Kümin
Mittelhäusern, 12 February 2018

Use of Sanosil Q-Jet Superior Together with Sanosil S015 to Fumigate Animal Rooms and Validation with Foot-and-Mouth Disease Virus (FMDV)

Dear Mr Rüegg

The Institute of Virology and Immunology (IVI) at Mittelhäusern, Switzerland, is the National Reference Laboratory for highly contagious animal diseases. The IVI diagnoses, researches and tests vaccines against a number of highly infectious pathogens, amongst those foot-and-mouth disease virus. The IVI is required to inactivate everything that leaves the facility. This may occur chemically (surface treatment or fumigation) or thermally. Also, animal rooms, especially, need to be fumigated following animal trials. The IVI decided to look for alternatives to formaldehyde fumigation and thus tested and validated the Q-Jet Superior with Sanosil S015. We were able to show that animal rooms can be efficiently, safely and successfully fumigated using the Q-Jet Superior. Besides using bacterial spore indicators to confirm successful fumigation, we also successfully tested and validated the fumigation with foot-and-mouth disease virus.

We can confirm that the Sanosil Q-Jet Superior with Sanosil S015 is a safe and efficient alternative to other fumigation methods on the market and can endorse it without reservations.

Yours sincerely,

Daniel Kümin, PhD
Study Lead

Kathrin Summermatter, PhD
Head of Biosafety

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