



HYGITEC Ltd. – Sofia, Bulgaria
15, Acad. Iv. Evst. Geshov, Blvd.
Fax/Phone: (0359 2) 954-11 19
www.hygitec.com :e-mail : info@hygitec.com



CERTIFICATE OF QUALITY
No. XAD2 -4 / 2017

SAMPLING TUBES FOR XAD-2 (2-HMP)

Characteristics	Value
Type of the sorbent tubes	Adsorbent Tubes XAD-2 (2-HMP)
Catalogue No.	ST 21-12
Lot #	18670
Expiration date	02.2019
Design parameters	
Overall length of the glass tube, mm	110 ± 2
Outside diameter, mm	6.0 ± 0.2
Internal diameter, mm	4.0 ± 0.2
Diameter of the tips, mm	2.0 ± 0.3
Weight of the back sorbent layer, mg	60 ± 1
Weight of the front sorbent layer, mg	120 ± 10
Characteristics of materials	
Glass tubes	Neutral glass
Type of glass wool separators	Silanized
Filling material	XAD-2 (2-HMP)

Production Manager: *.....*



Sofia, 2017

/Eng.-Chem. J. Vassileva/